

# **S12/S12W** PC ECG

With a data acquisition module and a set of comprehensive software program for acquiring, analyzing, storing and reviewing 12-channel ECGs, SS12 PC ECG system provides all-in-one solutions for both resting ECG and stress test.

# Wired and wireless solutions make your options more flexible

### Ability to work with both PC & Tablet

MEDIBLU'S PCG ECG offers you flexibility to view ECG traces on both PC and Tablet. (PC: Windows 7 32/64-bit, Windows Vista 32/64-bit. Windows XP SP2/SP3. Tablet: Windows 7 32/64-bit)

















Multi formats ECG reports export and can be sent by email to streamline your workflow.

### Wired Solution

Wired solution provides you with simplified installation.

#### Wireless Solution

Wireless solution provides the freedom of patients and clinicians during the stress test.





#### 12-channel Resting ECG

- Clear and reliable ECG waveforms acquisition
- Accurate measurements and interpretation testedwith CSE/MIT/AHA database

## Comprehensive Analysis

- Vector ECG, Timed VCG, Signal Average Electrocardiogram (SAECG)
- Heart Rate Turbulence (HRT), QT interval Dispersion
- High Frequency ECG, Frequency Spectrum cardiogram

#### Versatile Report

- PDF, FDA-XML, SCP, DICOM, JPG, BMP and DOC format reports transferrable via E-mail
- Compatible with GDT systems

#### Stress test

- Wired and wireless solutions
- Compatible with reputable treadmills and bicycles
- Support exercise blood pressure monitor
- Customized exercise protocol in addition to Bruce and Modified Bruce

PC-based ECG S12/S12W offers you an affordable and portable12-channel ECG solution with complete data analysis functions.











**QT Interval Dispersion** 

**HF ECG** 

**SAECG**